



ELK Biotechnology

RFP-Tag Mouse mAb

Catalog NO.: EM1010

For research use only.

Overview

Product name	RFP-Tag Mouse Monoclonal antibody
Source	Mouse
Applications	WB
Species reactivity	N/A
Recommended dilutions	WesternBlot:1/5000 NOTE: Optimal dilutions should be determined by the end user.
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	N/A
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	N/A
Background	The red fluorescent protein cloned from Discosoma coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein(GFP) from the Aequorea jellyfish.

